

INTRODUCING



Solutions for Neurophysiology



Wireless, Portable NIRS Imager



Advanced, Wireless Amplifier

The best portable and mobile, multimodal EEG-fNIRS setup

Cap
Compatible with both EEG electrodes and fNIRS optodes for seamless Integration.




NIRSport 2
Portable, wireless and truly whole-head, high density instrument featuring 23 to 200+ channels.

LIVEAMP
Wearable, wireless EEG amplifier that is expandable from 8 to 64 channels.



LiveAmp and NIRSport2 integrate with LSL for triggering and event synchronization for EEG fNIRS co-registration.

Solutions for EEG and fNIRS including a complete set of accessories is available through Symbiotic Devices.



1300 934 947

team@symbioticdevices.com.au

