

Arc Apollo[™] EEG

REMOTE EEG MONITORING SOLUTION

Arc Apollo ambulatory EEG allows physicians and technicians the ability to remotely view and monitor the quality of continuous ambulatory EEG data with HD Q-Video Mobile 3 with automatic infrared switching.



Arc Apollo[™] EEG

REMOTE EEG MONITORING SOLUTION

Apollo ambulatory EEG is designed for patient comfort and high-quality EEG. The recorder and amplifier weigh only 1 lb. 3 oz. and provide up to 96 hours of uninterrupted recording on a single charge, programmable patient event buttons, and continuous impedance checking.

Cadwell's remote EEG monitoring solution includes the Apollo Recorder, 32-channel Apollo Amplifier, Q-Video[®] Mobile 3 HD camera with automatic infrared switching, patient event microphone, a tablet or PC, cables, a case, and extra batteries. Add a mobile hotspot and remote access software for a complete remote monitoring solution.

Contact Symbiotic Devices today to schedule a demo.



Unit 6, 105-111 Ricketts Rd,
Mount Waverley, VIC 3149
ph 1300 934 947 | fx 1300 734 712
www.symbioticdevices.com.au
sales@symbioticdevices.com.au

