

Advanced Mobile EEG Solution



Outstanding data quality, 8-256 channels

Mobile and stationary equipment

Electrodes: Passive - Active - Active Dry

Intuitive, user-friendly software with direct access to data

Great match for various application fields

LiveAmp wireless amplifier for mobile EEG

- 8, 16, 32, 64 channel, 24-bit amplifier.
- Use an SD card or hard disk for data storage.
- Range of up to 30m
- Sensor extension available.



Up to 64 channels

actiCHamp (up to 160 channels)

High performance active 64 channel EEG amplifier. Active impedance feedback, 8 AUX inputs, 24-bit resolution, up to 160 channels, high sampling rate (up to 100KHz).

actiCap Active slim and snap electrodes

Generate excellent data under all recording conditions with live impedance feedback. Slim, lightweight and robust design.



1300 934 947

sales@symbioticdevices.com.au