


Sponge-based electrode system



- ✓ Ultra-fast measurement preparation
- ✓ Comfortable high-density EEG recordings
- ✓ Highly convenient - no gel residues
- ✓ Ag/AgCl sensors with exchangeable sponge tips
- ✓ Available with 32, 64, 96, 128 or 160 channels



Great match for time-restricted recordings

- ✓ If your research requires ultra-fast application, minimal invasiveness and maximum comfort, the R-Net is the perfect match.
- ✓ It is also the only sponge-based electrode system available that follows the international 10-10 positioning system.

New addition to our existing electrode types

- ✓ It expands our wide range of already available electrode cap types, i.e. passive gel, passive solid gel, active gel and active dry.

Fast preparation time

- ✓ It can be prepared with our potassium chloride electrolyte solution, allowing for the fastest preparation time out of all electrode systems.

Compatible with our mobile and stationary amplifiers

- ✓ It is compatible with our mobile LiveAmp systems, as well as our stationary gold standard BrainAmp amplifier series and our new versatile all-in-one lab amplifier, the actiCHamp Plus.

Quick and easy replacement of damaged electrodes

- ✓ Available repair kits allow you to fix a broken electrode yourself without the need to replace the whole R-Net.

Authorized Distributor



p 1300 934 947 **f** 1300 734 712
w www.symbioticdevices.com.au
e team@symbioticdevices.com.au
a Unit 6, 105-111 Ricketts Road
Mount Waverley, VIC 3149



www.brainproducts.com