



# IONM Electrodes, Probes and Accessories



## Single, Corkscrew, **Subdermal Needle Electrodes**

Subdermal Needle Electrodes are made of high tensile stainless steel and are ideal for both stimulating as well as recording. Symbiotic Devices offers needle electrodes in a great variety of lead wire length, styles and colours either in singles or pairs.

### Single, Disposable, Corkscrew, Stainless Steel, Needle Electrodes

Corkscrew subdermal needle electrodes are designed to stay securely in place while providing great signal.



SD Part Number	Description	Needle Ø	Lead Length	Qty
SDRL-NE1029	6-Colours Single	0.6mm / 23G	1.0m / 39"	12 / Box
SDRL-NE1030	6-Colours Single	0.6mm / 23G	1.5m / 59"	12 / Box
SDRL-NE1031	6-Colours Single	0.6mm / 23G	2.5m / 98"	12 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ

### Single, Disposable, Corkscrew, Subdermal Needle Electrodes - Generic

SD Part Number	Description	Needle Ø	Wire	Qty
SDFM-NE0054	6-Colours Single - Blue hub	0.6mm / 23G	1.2m / 47"	24 / Box
SDFM-NE0055	6-Colours Single - Blue hub	0.6mm / 23G	1.5m / 59"	24 / Box
SDFM-NE0056	6-Colours Single - Blue hub	0.6mm / 23G	2.0m / 79"	24 / Box

ARTG ID: 311383



## Twisted, Corkscrew, **Subdermal Needle Electrodes**

### Twisted Pair, Disposable, Corkscrew, Needle Electrodes - Generic

SD Part Number	Description	Needle Ø	Lead Length	Qty
SDFM-NE0057	Twisted Pair - 12 Colours	0.6mm / 23G	2.0m / 79"	24 / Box

ARTG ID: 311383



# Single, Subdermal Needle Electrodes

## Single, Disposable, Subdermal Needle Electrodes

SD Part Number	Description	Needle Length	Needle Ø	Lead Length	Qty
SDRL-NE0004	6-Colours Group 1	7mm	0.4mm	1.5m	24 / Box
SDRL-NE0005	6-Colours Group 1	7mm	0.4mm	2.5m	24 / Box
SDRL-NE1004	6-Colours Group 2	7mm	0.4mm	1.5m	24 / Box
SDRL-NE1005	6-Colours Group 2	7mm	0.4mm	2.5m	24 / Box
SDRL-NE0006	6-Colours Group 1	13mm	0.4mm	1.0m	24 / Box
SDRL-NE0007	6-Colours Group 1	13mm	0.4mm	1.5m	24 / Box
SDRL-NE0008	6-Colours Group 1	13mm	0.4mm	2.5m	24 / Box
SDRL-NE1009	6-Colours Group 2	13mm	0.4mm	2.5m	24 / Box
SDRL-NE0104	6-Colours Group 1	19mm	0.4mm	1.5m	24 / Box
SDRL-NE0012	6-Colours Group 1	19mm	0.4mm	2.5m	24 / Box
SDRL-NE0032	6-Colours Group 1	22mm	0.4mm	1.5m	24 / Box
SDRL-NE0033	6-Colours Group 1	22mm	0.4mm	2.5m	24 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ



SDRL-NE0004



7mm



13mm



19mm



22mm

## Single, Disposable, Subdermal Needle Electrodes - Generic

SD Part Number	Description	Needle Length	Needle Ø	Lead Length	Qty
SDFM-NE0031	6-Colours Single	7mm	0.35mm	1.5m	24 / Box
SDFM-NE1032	6-Colours Single	7mm	0.4mm	1.5m	24 / Box
SDFM-NE0035	6-Colours Single	13mm	0.4mm	1.0m	24 / Box
SDFM-NE0036	6-Colours Single	13mm	0.4mm	1.5m	24 / Box
SDFM-NE0037	6-Colours Single	13mm	0.4mm	2.0m	24 / Box
SDFM-NE0038	6-Colours Single	13mm	0.4mm	2.5m	24 / Box
SDFM-NE0032	Black wire only	15mm	0.35mm	1.0m	24 / Box
SDFM-NE0052	6-Colours Single	17mm	0.45mm	1.5m	24 / Box
SDFM-NE0051	6-Colours Single	23mm	0.4mm	1.5m	24 / Box

ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX



7mm Needle - 1.5m



13mm Needle - 1.5m

## Paired, Subdermal Needle Electrodes

### Paired, Disposable, Subdermal Needle Electrodes

SD Part Number	Length	Needle Ø	Lead Length	Colour	Qty
SDRL-NE0113	7mm	0.4mm	1.5m	6 - Group 1	10 / Box
SDRL-NE0114	7mm	0.4mm	2.5m	6 - Group 1	10 / Box
SDRL-NE0115	7mm	0.4mm	1.5m	6 - Group 2	10 / Box
SDRL-NE0116	7mm	0.4mm	2.5m	6 - Group 2	10 / Box
SDRL-NE0017	13mm	0.4mm	1.0m	6 - Group 1	10 / Box
SDRL-NE0018	13mm	0.4mm	1.5m	6 - Group 1	10 / Box
SDRL-NE0019	13mm	0.4mm	2.5m	6 - Group 1	10 / Box
SDRL-NE0020	13mm	0.4mm	1.0m	6 - Group 2	10 / Box
SDRL-NE0021	13mm	0.4mm	1.5m	6 - Group 2	10 / Box
SDRL-NE0022	13mm	0.4mm	2.5m	6 - Group 2	10 / Box
SDRL-NE0023	19mm	0.4mm	1.5m	6 - Group 1	10 / Box
SDRL-NE0024	19mm	0.4mm	2.5m	6 - Group 1	10 / Box
SDRL-NE0025	19mm	0.4mm	1.5m	6 - Group 2	10 / Box
SDRL-NE0026	19mm	0.4mm	2.5m	6 - Group 2	10 / Box
SDRL-NE0034	22mm	0.4mm	1.5m	6 - Group 1	10 / Box
SDRL-NE0035	22mm	0.4mm	2.5m	6 - Group 1	10 / Box
SDRL-NE0036	22mm	0.4mm	1.5m	6 - Group 2	10 / Box
SDRL-NE0037	22mm	0.4mm	2.5m	6 - Group 2	10 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ



13mm Paired



## Parallel, Subdermal Needle Electrodes

Parallel, Disposable, Subdermal Needle Paired Electrodes  
 Vacuum-smelted, high-elastic and medical grade stainless steel wires with Safe-touch DIN 42802



SD Part Number	Length	Needle Ø	Lead Length	Qty
SDC-NE0114	13mm	0.4mm / 30G	1.5m / 59"	10 / Box
SDC-NE0116	13mm	0.4mm / 30G	2.5m / 98"	10 / Box

ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX



Six-Parallel, Disposable, Subdermal Needle Electrodes

SD Part Number	Length	Needle Ø	Lead Length	Qty
SDC-NE0045	13mm	0.4mm / 30G	1.5m / 59"	10 / Box

ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX



## Twisted, Subdermal Needle Electrodes



### Twisted, Disposable, Subdermal Needle Electrodes

SD Part Number	Length	Needle Ø	Lead Length	Colour	Qty
SDRL-NE0028	13mm	0.4mm	2.5m	6 - Group 1	10 / Box
SDRL-NE0027	13mm	0.4mm	2.5m	6 - Group 2	10 / Box
SDRL-NE0038	13mm	0.4mm	2.5m	6 - Group 3	10 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ



### Twisted Pair, Disposable, Subdermal Needle Electrodes - Generic

SD Part Number	Length	Needle Ø	Lead Length	Colour	Qty
SDFM-NE0034	7mm	0.4mm	1.5m	12 Colours	12 Pair / Box
SDFM-NE0039	13mm	0.4mm	1.0m	12 Colours	12 Pair / Box
SDFM-NE0040	13mm	0.4mm	1.5m	12 Colours	12 Pair / Box
SDFM-NE0042	13mm	0.4mm	2.0m	12 Colours	12 Pair / Box
SDFM-NE0041	13mm	0.4mm	2.0m	12- Red & Black only	12 Pair / Box
SDFM-NE0043	13mm	0.4mm	2.5m	12 Colours	12 Pair / Box
SDFM-NE0044	15mm	0.4mm	1.5m	12 Colours	12 Pair / Box
SDFM-NE0139	15mm	0.4mm	2.0m	12 Colours	12 Pair / Box
SDFM-NE0140	17mm	0.4mm	2.0m	12 Colours	12 Pair / Box
SDFM-NE0141	17mm	0.4mm	1.5m	12 Colours	12 Pair / Box

ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX



13mm Twisted Pair



### 3-Twisted Set, Disposable, Subdermal Needle Electrodes - Generic



SD Part Number	Length	Needle Ø	Lead Length	Colour	Qty
SDFM-NE0030	20mm	0.3mm	1.0m	Red, Black, Green	12 / Box

ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX



### 5-Twisted Set, Disposable, Subdermal Needle Electrodes



SD Part Number	Needle Length	Needle Ø	Lead Length	Qty
SDRL-NE0127	13mm	0.4mm	1.5m	4 / Box
SDRL-NE0128	13mm	0.4mm	2.5m	4 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ



## Twisted, Subdermal Needle Electrodes

6-Twisted Set, Disposable, Subdermal Needle Electrodes 



SD Part Number	Needle Length	Needle Ø	Lead Length	Qty
SDC-NE0300	13mm	0.4mm / 27G	1.5m / 59"	10 / Box

ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX




6-Twisted Set, Disposable, Subdermal Needle Electrodes 



SD Part Number	Needle Length	Needle Ø	Lead Length	Qty
SDRL-NE0129	13mm	0.4mm	1.5m / 59"	3 / Box
SDRL-NE0130	13mm	0.4mm	2.5m / 98"	3 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ



## Other types, Subdermal Needle Electrodes

Bent, Disposable, Subdermal Needle Electrodes

A specially designed hook is ideal for long recordings, while bending angle improves skin contact and lower impedance of the electrode.



SD Part Number	Description	Angle	Needle	Lead	Qty
SDRL-NE0106	6-Colours	35°	13 x 0.4mm	1.0m / 39"	24 / Box
SDRL-NE0013	6-Colours	35°	13 x 0.4mm	2.5m / 98"	24 / Box
SDRL-NE0015	12-Colours	35°	13 x 0.4mm	1.5m / 59"	24 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ




  
6-Colours Group 1

ENT Dual Needle Electrodes, Disposable



SD Part Number	Needle Length	Needle Ø	Lead Length	Qty
SDFM-NE0019	12mm	0.4mm / 30G	1.5m / 59"	12 / Box

ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX



Guardian Subdermal Needle Electrode, Disposable

Protects the patient and medical personnel from accidental needle injury as well as delivering equivalent recording and stimulating quality.



SD Part Number	Description	Needle Length	Lead Length	Qty
SDRL-NE0900	Needle with Bandage in Protective Hub	13mm	1.5m	24 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ




## Stimulating Probes, Disposable

Designed for surgeons and specific applications. Made from stainless steel with handles.

### Bipolar Probes

SD Part Number	Description	Probe Length	Tip Ø	Qty
SDFM-P0002	Bipolar Direct Nerve Stimulator Probe	90mm	0.8mm	10 / Box

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA




Bipolar Direct

### Concentric Probes

SD Part Number	Description	Probe Length	Tip Ø	Qty
SDFM-P0007	Concentric Direct	90mm	1.3mm	10 / Box
SDFM-P0008	Bent Concentric Direct, Bent Tip Length 6mm	80 - 120mm	1.3mm	Each

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA



Concentric Direct Bent

### Flush Tip Stimulating Probes

SD Part Number	Description	Probe Length	Tip Ø	Lead	Qty
SDFM-P0003	Flush Tip Direct Nerve	80 - 150mm	0.8mm	2.0m	10 / Box
SDRL-P0006	FlushTip SlideShaft®	130mm - 270mm	0.2mm	2.0m	5 / Pk

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA

ARTG ID: 290956 | WAND: 170623-WAND-6OE2PQ



SlideShaft®  
Flush Tip

### Monopolar Stimulating Probes

SD Part Number	Description	Probe Length	Probe Ø	Tip Ø	Lead	Qty
SDRL-P0008	Extended	160mm	3.5mm	1.6mm	2.5m	Each
SDRL-P0011	With Removable Handle	200mm	1.5mm	1.5mm	2.5m	Each

ARTG ID: 290956 | WAND: 170623-WAND-6OE2PQ



Extended  
Monopolar

### Monopolar Stimulating Probes - Generic

SD Part Number	Description	Size	Qty
SDC-P0003	Probe, Monopolar Stimulating, Disposable	90mm	Each
SDC-P0004	Probe, Monopolar Stimulating, Disposable	200mm	10 / Box
SDC-P0005	Probe, Monopolar Stimulating, Disposable	275mm	5 / Box

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA



Monopolar  
Stimulating



# Stimulating Probes, Disposable

## Double Hook Stimulating Probes

SD Part Number	Probe Length	Probe Ø	Hook Degree	Shaft Spacing	Qty
SDC-P0016	90mm	1.2mm	90°	5mm	Each
SDFM-P0009	70-120mm	1.2mm	100°	5mm	Each
SDFM-P0010	70-120mm	1.2mm	180°	5mm	Each
SDFM-P0011	Right Angle 70-120mm	1.2mm	90°	5mm	Each

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA



Double Hook

## Triple Hook Stimulating Probes

SD Part Number	Probe Length	Probe Ø	Hook Degree	Shaft Spacing	Qty
SDFM-P0013	70-120mm	1.2mm	100°	5mm	Each
SDFM-P0014	70-120mm	1.2mm	180°	5mm	Each

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA



Triple Hook 100°




Triple Hook 180°

## Needle Stimulating Probes

SD Part Number	Needle Length	Needle Ø	Qty
SDFM-P0017	13mm	0.25mm	Each

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA




Needle Stimulator Probe

## Ball Tip Probes

### Ball Tip, Bent Direct Nerve Stimulator Probes

SD Part Number	Probe Length	Tip Ø	Lead	Qty
SDC-P0007	200mm	2.3mm	2m	Each

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA



Bent Nerve Stimulating Probe

### Ball Tip, Double Ball, Nerve Stimulating Probes

SD Part Number	Probe Length	Tip Ø	Lead	Qty
SDMI-P0015	80mm	2.0mm	2m	Each

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA




Double Ball Tip

### Ball Tip, Monopolar Probes

SD Part Number	Probe Length	Tip Ø	Lead	Qty
SDC-P0010	90 - 150mm	2.3mm	2m	10 / Box
SDC-P0011	160mm / 6"	2.3mm	2m	10 / Box
SDC-P0012	200mm / 8"	2.3mm	2m	10 / Box

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA



Ball Tip Monopolar

### Ball Tip, SlideShaft® Stimulating Probes

A neuro stimulating probe with adjustable shaft length according to surgical procedure.



SD Part Number	Probe Length	Probe Ø	Tip Ø	Lead	Qty
SDRL-P0007	270mm Max.	1.85 - 2.5mm	2.3mm	2.5m	Each

ARTG ID: 290956 | WAND: 170623-WAND-6OE2PQ



SlideShaft® Ball Tip Probe

### Ball Tip, Standard Stimulating Probes



SD Part Number	Probe Length	Probe Ø	Tip Ø	Lead	Qty
SDRL-P0001	100mm	1.5mm	3mm	2.5m	Each
SDRL-P0009	100mm	1.5mm	3mm	2.5m	10 / Box

ARTG ID: 290956 | WAND: 170623-WAND-6OE2PQ




Standard Ball Tip

### Ball Tip, Tapered Stimulating Probes



SD Part Number	Probe Length	Probe Ø	Tip Ø	Lead	Qty
SDRL-P0003	150mm	1.85 - 2.5mm	2.3mm	2.5m	Each
SDRL-P0002	150mm	1.85 - 2.5mm	2.3mm	2.5m	10 / Box
SDRL-P0004	200mm	1.85 - 2.5mm	2.3mm	2.5m	10 / Box
SDRL-P0005	200mm	1.85 - 2.5mm	2.3mm	2.5m	Each

ARTG ID: 290956 | WAND: 170623-WAND-6OE2PQ



Tapered Ball Tip

## Accessories

### Clips

SD Part Number	Description	Qty
SDC-CC0032	Clip-On Lead wire, Plunger Style	10 / Box
SDC-CC0033	Multi-Stage Clip-on Lead wire	10 / Box
ARTG ID: 356742		
SDRL-MI0001	Alligator Clip (Sterile), Disposable, Lead Length 2.5m	12 / Box

ARTG ID: 274257 | WAND: 170623-WAND-6OE2PQ



Multi-Stage Clip



## Surface Electrodes, Disposable

Surface Electrodes, Single Pre-gelled, with Lead wires, Disposable  
Four mil carbon film, silver/silver chloride coated.

SD Part Number	Wire Colours	Electrode Size	Lead Length	Qty
SDC-SE0050	Red, Blue, Black, Green	27 x 20mm	2.5m	12 / Box
SDFM-SE0039	White	15 x 20mm	1.0m / 39"	144 / Box

ARTG ID: 310991 | WAND: 181026-WAND-6RB3UX



SDFM-SE0039

Surface Electrodes, Twisted Pre-gelled, with Lead wires, Disposable  
Four mil carbon film, silver/silver chloride coated.

SD Part Number	Description	Electrode Size	Lead Length	Qty
SDC-SE0051	Twisted Pre-gelled	20 x 15mm	1.5m	24 / Box
SDC-SE0052	Twisted Pre-gelled	20 x 15mm	2.5m	24 / Box
SDC-SE0053	Twisted Pre-gelled	27 x 20mm	1.5m	24 / Box
SDC-SE0054	Twisted Pre-gelled	27 x 20mm	2.5m	24 / Box

ARTG ID: 310991 | WAND: 181026-WAND-6RB3UX




## Surface Electrode, Ground, with Lead wire, Disposable

SD Part Number	Description	Electrode Size	Lead Length	Qty
SDC-SE1000	Ground, Pre-gelled	50 x 40mm	1.0m / 39"	20 / Box
SDC-SE1001	Ground, Pre-gelled	50 x 40mm	1.5m / 59"	20 / Box
SDC-SE1002	Ground, Pre-gelled	50 x 40mm	2.0m / 79"	20 / Box
SDC-SE1003	Ground, Pre-gelled	50 x 40mm	2.5m / 98"	20 / Box
SDC-SE1004	Ground, Pre-gelled	50 x 40mm	3.0m / 118"	20 / Box

ARTG ID: 310991 | WAND: 181026-WAND-6RB3UX

See next page for other sizes >




## Surface Electrodes, Disposable



### Surface Electrodes, Ground Sticky Pad, Disposable

SD Part Number	Description	Size	Lead Length	Qty
SDRL-SE0045	Single, Pre-gelled, Green	45 x 35mm	1.0m	40 / Box
SDRL-SE0046	Single, Pre-gelled, Green	45 x 35mm	1.5m	40 / Box
SDRL-SE0047	Single, Pre-gelled, Green	45 x 35mm	2.5m	40 / Box

ARTG ID: 290956 | WAND: I70623-WAND-6OE2PQ



### Surface Electrodes, Single Stimulating Sticky Pads, Disposable

SD Part Number	Description	Size	Lead Length	Qty
SDRL-SE0042	Small, Pre-gelled, Red, Blue, Black, Green	20 x 15mm	1.5m	80 / Box
SDRL-SE0043	Small, Pre-gelled, Red, Blue, Black, Green	20 x 15mm	2.5m	80 / Box
SDRL-SE0044	Small, Pre-gelled, Red, Blue, Black, Green	20 x 15mm	3.0m	80 / Box
SDRL-SE0040	Large, Pre-gelled, Red, Blue, Black, Green	20 x 25mm	1.5m	80 / Box
SDRL-SE0049	Large, Pre-gelled, Red, Blue, Black, Green	20 x 25mm	2.5m	80 / Box

ARTG ID: 290957 | WAND: I70623-WAND-6OE2PQ



### Surface Electrodes, Parallel Pair Stimulating Sticky Pads, Disposable

SD Part Number	Description	Size	Lead Length	Qty
SDRL-SE0060	Small, Pre-gelled, Black / Red, Blue / Green	20 x 15mm	3.0m	20 Pair / Box
SDRL-SE0061	Small, Pre-gelled, Black / Red, Blue / Green	20 x 15mm	1.5m	20 Pair / Box
SDRL-SE0062	Small, Pre-gelled, Black / Red, Blue / Green	20 x 15mm	2.5m	40 Pair / Box
SDRL-SE0063	Large, Pre-gelled, Black / Red, Blue / Green	20 x 25mm	1.5m	40 Pair / Box
SDRL-SE1064	Large, Pre-gelled, Black / Red, Blue / Green	20 x 25mm	2.5m	40 Pair / Box

ARTG ID: 290957 | WAND: I70623-WAND-6OE2PQ



# Surface Electrodes, Disposable



## Surface Electrodes, Recording Sticky Pads, Single, Disposable

SD Part Number	Description	Size	Lead Length	Qty
SDRL-SE0050	Small, Pre-gelled, Red, Blue, Black, Green	20 x 15mm	1.0m	80 / Box
SDRL-SE0051	Small, Pre-gelled, Red, Blue, Black, Green	20 x 15mm	1.5m	80 / Box
SDRL-SE0052	Small, Pre-gelled, Red, Blue, Black, Green	20 x 15mm	2.5m	80 / Box
SDRL-SE0053	Large, Pre-gelled, Red, Blue, Black, Green	20 x 25mm	1.5m	80 / Box
SDRL-SE0054	Large, Pre-gelled, Red, Blue, Black, Green	20 x 25mm	2.5m	80 / Box

ARTG ID: 290957 | WAND: 170623-WAND-6OE2PQ

## Surface Electrodes, Recording Sticky Pads, Parallel Pair, Disposable

SD Part Number	Description	Size	Lead Length	Qty
SDRL-SE0056	Small, Pre-gelled, Black / Red, Blue / Green	20 x 15mm	1.5m	20 Pair / Box
SDRL-SE0057	Small, Pre-gelled, Black / Red, Blue / Green	20 x 15mm	2.5m	20 Pair / Box
SDRL-SE0058	Small, Pre-gelled, Black / Red, Blue / Green	20 x 25mm	1.5m	40 Pair / Box
SDRL-SE0059	Large, Pre-gelled, Black / Red, Blue / Green	20 x 25mm	2.5m	40 Pair / Box

ARTG ID: 290957 | WAND: 170623-WAND-6OE2PQ



# Extension Cables

SD Part Number	Description	Length	Qty
SDC-CC0023	Cadwell Dual Extension Lead Wire, 60" ● ○	152cm	Each
SDC-CC0024	Cadwell Dual Extension Lead Wire, 180" ● ○	457cm	Each
SDC-CC0025	Cadwell Quad Extension Lead Wire, 72" ● ● ● ●	183cm	Each
SDC-CC0030	Reference Jumper. Two stackable 1.5mm, ultra-flexible ●	15cm	Each
SDC-CC0029	Reference Jumper. Two stackable 1.5mm, ultra-flexible ●	249cm	Each
ARTG ID: 356742			
SDFM-CC0019	Extension Lead Wire 6 Colours ● ● ● ● ● ●	100cm	6 / Pk
SDFM-CC0021	Extension Lead Wire 6 Colours ● ● ● ● ● ●	150cm	6 / Pk
SDFM-CC0023	Extension Lead Wire 6 Colours ● ● ● ● ● ●	200cm	6 / Pk

ARTG ID: 310991



Dual Extension Wire



SDC-CC0030



Generic Extension Lead Wire

\* ARTG ID: 356742  
\*\* ARTG ID: 310991



Symbiotic Devices was founded in 2012 on a passion not only for neuroscience and neurotechnology, but also to provide clients with uncompromised quality, reliability and assistance. Striving to contribute to a sustainable successful research and clinical sectors across Australia and New Zealand, Symbiotic Devices aims to balance quality, accessibility and customisation.

We know that each neurophysiology project will present unique challenges, and our wide range of partners means we are able to provide customised and dynamic solutions to fit our client's needs.

First partnering with Brain Products, EasyCap and Rogue Research, Symbiotic Devices now proudly represents CED, Deymed, Cadwell, NeuroConn/NeuroCare, Rhythmlink, Pearl Technology and VPixx.

If you want the best solutions for your neuroscience and neurophysiology needs, you can trust Symbiotic Devices to provide the highest quality equipment, applications and support while contributing to a sustainable and ethical future for neuroscience practices.



**SYMBIOTIC**  
DEVICES

**au 1300 934 947 nz 09 886 4877**  
[team@symbioticdevices.com.au](mailto:team@symbioticdevices.com.au)

Unit 6, 105-111 Ricketts Rd,  
Mount Waverley, VIC 3149

[www.symbioticdevices.com.au](http://www.symbioticdevices.com.au)



**Contact us**

[team@symbioticdevices.com.au](mailto:team@symbioticdevices.com.au) | Ph: 1300 934 947