



Australian Government
Department of Health
Therapeutic Goods Administration

Public Summary

Summary for ARTG Entry: 311382 Symbiotic Devices Pty Ltd - Nerve-locating probe, single-use

ARTG entry for Medical Device Included Class IIa
Sponsor Symbiotic Devices Pty Ltd
Postal Address Suite 138 / 245 Thomas Street, DANDENONG, VIC, 3175
Australia
ARTG Start Date 10/11/2018
Product Category Medical Device Class IIa
Status Active
Approval Area Medical Devices

Conditions

- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods Act 1989 and Part 5, Division 5.2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information.
- Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.

Manufacturers

Name	Address
Xi' an Friendship Medical Electronics Co Ltd	No9 Gao Xin 1st Road High-Tech Development Zone Xi an, Shaanxi, 710075 China

Products

1 . Nerve-locating probe, single-use

Product Type	Single Device Product	Effective Date	10/11/2018
GMDN	61089 Nerve-locating probe, single-use		
Intended Purpose	A sterile, hand-held surgical instrument intended to function as an electrical conductor to locate a nerve during open surgery, and for diagnostic purposes (EMG), by delivering controlled electrical impulses to a specific body site for the measurement of a nerve response.		

Specific Conditions

No Specific Conditions included on Record

© Commonwealth of Australia. This work is copyright. You are not permitted to re-transmit, distribute or commercialise the material without obtaining prior written approval from the Commonwealth. Further details can be found at <http://www.tga.gov.au/about/website-copyright.htm>.

Public Summary