

IONM Electrodes, Probes and Accessories - NZ



Single, Corkscrew, Subdermal Needle Electrodes

Subdermal Needle Electrodes are made of high tensile stainless steel and are ideal for both stimulating as well as recording. Symbiotic Devices offers needle electrodes in a great variety of lead wire length, styles and colours either in singles or pairs.

Single, Disposable, Corkscrew, Stainless Steel, Needle Electrodes

Corkscrew subdermal needle electrodes are designed to stay securely in place while providing great signal.

SD Part Number	Description	Needle Ø	Lead Length	Qty
SDRL-NE1028	4-Colours Single	0.6mm / 23G	1.5m / 59"	12 / Box
SDRL-NE1029	6-Colours Single	0.6mm / 23G	1.0m / 39"	12 / Box
SDRL-NE1030	6-Colours Single	0.6mm / 23G	1.5m / 59"	12 / Box
SDRL-NE1031	6-Colours Single	0.6mm / 23G	2.5m / 98"	12 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ



Single, Disposable, Corkscrew, Subdermal Needle Electrodes - Generic

SD Part Number	Description	Needle Ø	Wire	Qty
SDFM-NE0054	6-Colours Single - Blue hub	0.6mm / 23G	1.2m / 47"	24 / Box
SDFM-NE0055	6-Colours Single - Blue hub	0.6mm / 23G	1.5m / 59"	24 / Box
SDFM-NE0056	6-Colours Single - Blue hub	0.6mm / 23G	2.0m / 79"	24 / Box



ARTG ID: 311383



Twisted, Corkscrew, Subdermal Needle Electrodes

Twisted, Disposable, Corkscrew, Stainless Steel, Needle Electrodes

Corkscrew twisted needle electrodes offer the quality of single needle electrodes but through twisted colour-coded lead wires.



SD Part Number	Description	Needle Ø	Lead Length	Qty
SDRL-NE1026	2-Colours Twisted	0.6mm / 23G	1.5m / 59"	6 / Box
SDRL-NE1027	4-Colours Twisted	0.6mm / 23G	1.5m / 59"	3 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ

Twisted Pair, Disposable, Corkscrew, Needle Electrodes



SD Part Number	Description	Needle Ø	Lead Length	Qty
SDC-NE0155	Twisted Pair	0.6mm / 23G	1.2m / 47"	12 / Box
SDC-NE0156	Twisted Pair	0.6mm / 23G	1.5m / 59"	12 / Box
SDC-NE0157	Twisted Pair	0.6mm / 23G	2.0m / 79"	12 / Box
SDC-NE0158	Twisted Pair	0.6mm / 23G	2.5m / 98"	12 / Box

ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX

Twisted Pair, Disposable, Corkscrew, Needle Electrodes - Generic



SD Part Number	Description	Needle Ø	Lead Length	Qty
SDFM-NE0057	Twisted Pair - 12 Colours	0.6mm / 23G	2.0m / 79"	24 / Box

ARTG ID: 311383



Single, Subdermal Needle Electrodes

Single, Disposable, Subdermal Needle Electrodes



SDRL-NE0004



7mm



13mm



19mm



22mm

SD Part Number	Description	Needle Length	Needle Ø	Lead Length	Qty
SDRL-NE0004	6-Colours Group 1	7mm	0.4mm	1.5m	24 / Box
SDRL-NE0005	6-Colours Group 1	7mm	0.4mm	2.5m	24 / Box
SDRL-NE1004	6-Colours Group 2	7mm	0.4mm	1.5m	24 / Box
SDRL-NE1005	6-Colours Group 2	7mm	0.4mm	2.5m	24 / Box
SDRL-NE0006	6-Colours Group 1	13mm	0.4mm	1.0m	24 / Box
SDRL-NE0007	6-Colours Group 1	13mm	0.4mm	1.5m	24 / Box
SDRL-NE0008	6-Colours Group 1	13mm	0.4mm	2.5m	24 / Box
SDRL-NE1009	6-Colours Group 2	13mm	0.4mm	2.5m	24 / Box
SDRL-NE0104	6-Colours Group 1	19mm	0.4mm	1.5m	24 / Box
SDRL-NE0012	6-Colours Group 1	19mm	0.4mm	2.5m	24 / Box
SDRL-NE0032	6-Colours Group 1	22mm	0.4mm	1.5m	24 / Box
SDRL-NE0033	6-Colours Group 1	22mm	0.4mm	2.5m	24 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ

Single, Disposable, Subdermal Needle Electrodes - Generic

SD Part Number	Description	Needle Length	Needle Ø	Lead Length	Qty
SDFM-NE0031	6-Colours Single	7mm	0.35mm	1.5m	24 / Box
SDFM-NE1032	6-Colours Single	7mm	0.4mm	1.5m	24 / Box
SDFM-NE0035	6-Colours Single	13mm	0.4mm	1.0m	24 / Box
SDFM-NE0036	6-Colours Single	13mm	0.4mm	1.5m	24 / Box
SDFM-NE0037	6-Colours Single	13mm	0.4mm	2.0m	24 / Box
SDFM-NE0038	6-Colours Single	13mm	0.4mm	2.5m	24 / Box
SDFM-NE0032	Black wire only	15mm	0.35mm	1.0m	24 / Box
SDFM-NE0052	6-Colours Single	17mm	0.45mm	1.5m	24 / Box
SDFM-NE0051	6-Colours Single	23mm	0.4mm	1.5m	24 / Box

ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX



7mm Needle - 1.5m



13mm Needle - 1.5m

Paired, Subdermal Needle Electrodes

Paired, Disposable, Subdermal Needle Electrodes

SD Part Number	Length	Needle Ø	Lead Length	Colour	Qty
SDRL-NE0113	7mm	0.4mm	1.5m	6 - Group 1	10 / Box
SDRL-NE0114	7mm	0.4mm	2.5m	6 - Group 1	10 / Box
SDRL-NE0115	7mm	0.4mm	1.5m	6 - Group 2	10 / Box
SDRL-NE0116	7mm	0.4mm	2.5m	6 - Group 2	10 / Box
SDRL-NE0017	13mm	0.4mm	1.0m	6 - Group 1	10 / Box
SDRL-NE0018	13mm	0.4mm	1.5m	6 - Group 1	10 / Box
SDRL-NE0019	13mm	0.4mm	2.5m	6 - Group 1	10 / Box
SDRL-NE0020	13mm	0.4mm	1.0m	6 - Group 2	10 / Box
SDRL-NE0021	13mm	0.4mm	1.5m	6 - Group 2	10 / Box
SDRL-NE0022	13mm	0.4mm	2.5m	6 - Group 2	10 / Box
SDRL-NE0023	19mm	0.4mm	1.5m	6 - Group 1	10 / Box
SDRL-NE0024	19mm	0.4mm	2.5m	6 - Group 1	10 / Box
SDRL-NE0025	19mm	0.4mm	1.5m	6 - Group 2	10 / Box
SDRL-NE0026	19mm	0.4mm	2.5m	6 - Group 2	10 / Box
SDRL-NE0034	22mm	0.4mm	1.5m	6 - Group 1	10 / Box
SDRL-NE0035	22mm	0.4mm	2.5m	6 - Group 1	10 / Box
SDRL-NE0036	22mm	0.4mm	1.5m	6 - Group 2	10 / Box
SDRL-NE0037	22mm	0.4mm	2.5m	6 - Group 2	10 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ



13mm Paired



Parallel, Subdermal Needle Electrodes

Parallel, Disposable, Subdermal Needle Paired Electrodes

Vacuum-smelted, high-elastic and medical grade stainless steel wires with Safe-touch DIN 42802



SD Part Number	Length	Needle Ø	Lead Length	Qty
SDC-NE0113	13mm	0.4mm / 30G	1.0m / 39"	10 / Box
SDC-NE0114	13mm	0.4mm / 30G	1.5m / 59"	10 / Box
SDC-NE0115	13mm	0.4mm / 30G	2.0m / 79"	10 / Box
SDC-NE0116	13mm	0.4mm / 30G	2.5m / 98"	10 / Box

ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX



Six-Parallel, Disposable, Subdermal Needle Electrodes



SD Part Number	Length	Needle Ø	Lead Length	Qty
SDC-NE0045	13mm	0.4mm / 30G	1.5m / 59"	10 / Box

ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX



Twisted, Subdermal Needle Electrodes



Twisted, Disposable, Subdermal Needle Electrodes

SD Part Number	Length	Needle Ø	Lead Length	Colour	Qty
SDRL-NE0038	13mm	0.4mm	1.5m	6 - Group 3	12 / Box
SDRL-NE0028	13mm	0.4mm	2.5m	6 - Group 1	10 / Box
SDRL-NE0027	13mm	0.4mm	2.5m	6 - Group 2	10 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ



IONM ELECTRODES

Twisted Pair, Disposable, Subdermal Needle Electrodes - Generic

SD Part Number	Length	Needle Ø	Lead Length	Colour	Qty
SDFM-NE0034	7mm	0.4mm	1.5m	12 Colours	12 Pair / Box
SDFM-NE0039	13mm	0.4mm	1.0m	12 Colours	12 Pair / Box
SDFM-NE0040	13mm	0.4mm	1.5m	12 Colours	12 Pair / Box
SDFM-NE0042	13mm	0.4mm	2.0m	12 Colours	12 Pair / Box
SDFM-NE0041	13mm	0.4mm	2.0m	12- Red & Black only	12 Pair / Box
SDFM-NE0043	13mm	0.4mm	2.5m	12 Colours	12 Pair / Box
SDFM-NE0044	15mm	0.4mm	1.5m	12 Colours	12 Pair / Box
SDFM-NE0139	15mm	0.4mm	2.0m	12 Colours	12 Pair / Box
SDFM-NE0140	17mm	0.4mm	2.0m	12 Colours	12 Pair / Box
SDFM-NE0141	17mm	0.4mm	1.5m	12 Colours	12 Pair / Box
SDFM-NE0033	23mm	0.35mm	2.0m	12 Colours	12 Pair / Box

ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX



13mm Twisted Pair



3-Twisted Set, Disposable, Subdermal Needle Electrodes - Generic ● ● ●

SD Part Number	Length	Needle Ø	Lead Length	Colour	Qty
SDFM-NE0030	20mm	0.3mm	1.0m	Red, Black, Green	12 / Box


ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX



5-Twisted Set, Disposable, Subdermal Needle Electrodes ● ● ● ● ●

SD Part Number	Needle Length	Needle Ø	Lead Length	Qty
SDRL-NE0127	13mm	0.4mm	1.5m	4 / Box
SDRL-NE0128	13mm	0.4mm	2.5m	4 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ



Twisted, Subdermal Needle Electrodes

6-Twisted Set, Disposable, Subdermal Needle Electrodes 

 CADWELL®

SD Part Number	Needle Length	Needle Ø	Lead Length	Qty
SDC-NE0300	13mm	0.4mm / 27G	1.5m / 59"	10 / Box

ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX



6-Twisted Set, Disposable, Subdermal Needle Electrodes 

 RhythmLink®

SD Part Number	Needle Length	Needle Ø	Lead Length	Qty
SDRL-NE0129	13mm	0.4mm	1.5m / 59"	3 / Box
SDRL-NE0130	13mm	0.4mm	2.5m / 98"	3 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ



Other types, Subdermal Needle Electrodes

Bent, Disposable, Subdermal Needle Electrodes

A specially designed hook is ideal for long recordings, while bending angle improves skin contact and lower impedance of the electrode.

 RhythmLink®

SD Part Number	Description	Angle	Needle	Lead	Qty
SDRL-NE0106	6-Colours	35°	13 x 0.4mm	1.0m / 39"	24 / Box
SDRL-NE0013	6-Colours	35°	13 x 0.4mm	2.5m / 98"	24 / Box
SDRL-NE0015	12-Colours	35°	13 x 0.4mm	1.5m / 59"	24 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ




ENT Dual Needle Electrodes, Disposable

 CADWELL®

SD Part Number	Needle Length	Needle Ø	Lead Length	Qty
SDFM-NE0018	12mm	0.4mm / 30G	1.0m / 39"	12 / Box
SDFM-NE0019	12mm	0.4mm / 30G	1.5m / 59"	12 / Box
SDFM-NE0020	12mm	0.4mm / 30G	2.0m / 79"	12 / Box
SDFM-NE0021	12mm	0.4mm / 30G	2.5m / 98"	12 / Box

ARTG ID: 311383 | WAND: 181026-WAND-6RB3UX



Guardian Subdermal Needle Electrode, Disposable

Protects the patient and medical personnel from accidental needle injury as well as delivering equivalent recording and stimulating quality.

 RhythmLink®

SD Part Number	Description	Needle Length	Lead Length	Qty
SDRL-NE0900	Needle with Bandage in Protective Hub	13mm	1.5m	24 / Box

ARTG ID: 274256 | WAND: 170623-WAND-6OE2PQ




Stimulating Probes, Disposable

Designed for surgeons and specific applications. Made from stainless steel with handles.

Bipolar Probes

SD Part Number	Description	Probe Length	Tip Ø	Qty
SDFM-P0002	Bipolar Direct Nerve Stimulator Probe	90mm	0.8mm	10 / Box

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA




Bipolar Direct

Concentric Probes

SD Part Number	Description	Probe Length	Tip Ø	Qty
SDFM-P0007	Concentric Direct	90mm	1.3mm	10 / Box
SDFM-P0008	Bent Concentric Direct, Bent Tip Length 6mm	80 - 120mm	1.3mm	Each

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA




Concentric Direct Bent

Flush Tip Stimulating Probes

SD Part Number	Description	Probe Length	Tip Ø	Lead	Qty
SDFM-P0003	Flush Tip Direct Nerve	80 - 150mm	0.8mm	2.0m	10 / Box
SDRL-P0006	FlushTip SlideShaft®	130mm - 270mm	0.2mm	2.0m	5 / Pk

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA

ARTG ID: 290956 | WAND: 170623-WAND-6OE2PQ



SlideShaft®
Flush Tip

Monopolar Stimulating Probes



SD Part Number	Description	Probe Length	Probe Ø	Tip Ø	Lead	Qty
SDRL-P0008	Extended	160mm	3.5mm	1.6mm	2.5m	Each
SDRL-P0011	With Removable Handle	200mm	1.5mm	1.5mm	2.5m	Each

ARTG ID: 290956 | WAND: 170623-WAND-6OE2PQ



Extended
Monopolar

Monopolar Stimulating Probes - Generic

SD Part Number	Description	Size	Qty
SDC-P0003	Probe, Monopolar Stimulating, Disposable	90mm	Each
SDC-P0004	Probe, Monopolar Stimulating, Disposable	200mm	10 / Box
SDC-P0005	Probe, Monopolar Stimulating, Disposable	275mm	5 / Box
SDC-P0008	Probe, Monopolar Stimulating, Disposable	335mm	Each

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA



Monopolar
Stimulating

Stimulating Probes, Disposable

Double Hook Stimulating Probes

SD Part Number	Probe Length	Probe Ø	Hook Degree	Shaft Spacing	Qty
SDC-P0016	90mm	1.2mm	90°	5mm	Each
SDFM-P0009	70-120mm	1.2mm	100°	5mm	Each
SDFM-P0010	70-120mm	1.2mm	180°	5mm	Each
SDFM-P0011	Right Angle 70-120mm	1.2mm	90°	5mm	Each

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA



Double Hook

Triple Hook Stimulating Probes

SD Part Number	Probe Length	Probe Ø	Hook Degree	Shaft Spacing	Qty
SDFM-P0013	70-120mm	1.2mm	100°	5mm	Each
SDFM-P0014	70-120mm	1.2mm	180°	5mm	Each

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA



Triple Hook 100°



Triple Hook 180°

Needle Stimulating Probes

SD Part Number	Needle Length	Needle Ø	Qty
SDFM-P0015	13mm	0.4mm	Each
SDFM-P0016	13mm	0.3mm	Each
SDFM-P0017	13mm	0.25mm	Each

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA



Needle Stimulator Probe




Ball Tip Probes

Ball Tip, Bent Direct Nerve Stimulator Probes

SD Part Number	Probe Length	Tip Ø	Lead	Qty
SDC-P0007	200mm	2.3mm	2m	Each

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA



Bent Nerve Stimulating Probe

Ball Tip, Double Ball, Nerve Stimulating Probes

SD Part Number	Probe Length	Tip Ø	Lead	Qty
SDMI-P0015	80mm	2.0mm	2m	Each

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA



Double Ball Tip

Ball Tip, Monopolar Probes

SD Part Number	Probe Length	Tip Ø	Lead	Qty
SDC-P0010	90 - 150mm	2.3mm	2m	10 / Box
SDC-P0011	160mm / 6"	2.3mm	2m	10 / Box
SDC-P0012	200mm / 8"	2.3mm	2m	10 / Box

ARTG ID: 311382 | WAND: 181026-WAND-6RB3UA



Ball Tip Monopolar

Ball Tip, SlideShaft® Stimulating Probes

A neuro stimulating probe with adjustable shaft length according to surgical procedure.



SD Part Number	Probe Length	Probe Ø	Tip Ø	Lead	Qty
SDRL-P0007	270mm Max.	1.85 - 2.5mm	2.3mm	2.5m	Each

ARTG ID: 290956 | WAND: 170623-WAND-6OE2PQ



SlideShaft® Ball Tip Probe

Ball Tip, Standard Stimulating Probes



SD Part Number	Probe Length	Probe Ø	Tip Ø	Lead	Qty
SDRL-P0001	100mm	1.5mm	3mm	2.5m	Each
SDRL-P0009	100mm	1.5mm	3mm	2.5m	10 / Box

ARTG ID: 290956 | WAND: 170623-WAND-6OE2PQ



Standard Ball Tip

Ball Tip, Tapered Stimulating Probes



SD Part Number	Probe Length	Probe Ø	Tip Ø	Lead	Qty
SDRL-P0003	150mm	1.85 - 2.5mm	2.3mm	2.5m	Each
SDRL-P0002	150mm	1.85 - 2.5mm	2.3mm	2.5m	10 / Box
SDRL-P0004	200mm	1.85 - 2.5mm	2.3mm	2.5m	10 / Box
SDRL-P0005	200mm	1.85 - 2.5mm	2.3mm	2.5m	Each

ARTG ID: 290956 | WAND: 170623-WAND-6OE2PQ



Tapered Ball Tip

ET Tube Electrodes, Disposable



Single Channel Laryngeal Surface

SD Part Number	ET Tube Size	Lead Length	Qty
SDC-MI0003	6.0 - 7.0	80cm	Each
SDC-MI0004	7.5 - 10.0	80cm	Each



SDC-MI0004

Accessories

Clips

SD Part Number	Description	Qty
SDC-CC0032	Clip-On Lead wire, Plunger Style	10 / Box
SDC-CC0033	Multi-Stage Clip-on Lead wire	10 / Box
ARTG ID: 356742		
SDRL-MI0001	Alligator Clip (Sterile), Disposable, Lead Length 2.5m	12 / Box

ARTG ID: 274257 | WAND: 170623-WAND-6OE2PQ



Multi-Stage Clip



Disc Electrodes, Disposable




Disc Electrodes with Leads, Disposable

Four mil carbon film, silver/silver chloride coated, improved hydrogel.

SD Part Number	Wire Colours	Electrode Size	Lead Length	Qty
SDC-SE0043	White, Green, Black, Red	20mm / 3/4"	0.5m	48 / Box
SDC-SE0041	Green, Blue, Red, Black	20mm / 3/4"	1m	20 / Box
SDC-SE0044	White, Green, Black, Red	20mm / 3/4"	1m	48 / Box
SDC-SE0045	White, Green, Black, Red	20mm / 3/4"	2m	48 / Box
SDC-SE0046	White, Green, Black, Red	20mm / 3/4"	3m	48 / Box

ARTG ID: 310991 | WAND: 181026-WAND-6RB3UX




Surface Electrodes, Disposable

Surface Electrodes, Single Pre-gelled, with Lead wires, Disposable
Four mil carbon film, silver/silver chloride coated.



SD Part Number	Wire Colours	Electrode Size	Lead Length	Qty
SDC-SE0047	Red, Blue, Black, Green	20 x 15mm	1.5m	12 / Box
SDC-SE0048	Red, Blue, Black, Green	20 x 15mm	2.5m	12 / Box
SDC-SE0049	Red, Blue, Black, Green	27 x 20mm	1.5m	12 / Box
SDC-SE0050	Red, Blue, Black, Green	27 x 20mm	2.5m	12 / Box

ARTG ID: 310991 | WAND: 181026-WAND-6RB3UX



Surface Electrodes, Twisted Pre-gelled, with Lead wires, Disposable
Four mil carbon film, silver/silver chloride coated.

SD Part Number	Description	Electrode Size	Lead Length	Qty
SDC-SE0051	Twisted Pre-gelled	20 x 15mm	1.5m	24 / Box
SDC-SE0052	Twisted Pre-gelled	20 x 15mm	2.5m	24 / Box
SDC-SE0053	Twisted Pre-gelled	27 x 20mm	1.5m	24 / Box
SDC-SE0054	Twisted Pre-gelled	27 x 20mm	2.5m	24 / Box

ARTG ID: 310991 | WAND: 181026-WAND-6RB3UX




Surface Electrode, Ground, with Lead wire, Disposable

SD Part Number	Description	Electrode Size	Lead Length	Qty
SDC-SE1000	Ground, Pre-gelled	50 x 40mm	1.0m / 39"	20 / Box
SDC-SE1001	Ground, Pre-gelled	50 x 40mm	1.5m / 59"	20 / Box
SDC-SE1002	Ground, Pre-gelled	50 x 40mm	2.0m / 79"	20 / Box
SDC-SE1003	Ground, Pre-gelled	50 x 40mm	2.5m / 98"	20 / Box
SDC-SE1004	Ground, Pre-gelled	50 x 40mm	3.0m / 118"	20 / Box

ARTG ID: 310991 | WAND: 181026-WAND-6RB3UX

See next page for other sizes >




Surface Electrodes, Disposable



Surface Electrodes, Ground Sticky Pad, Disposable

SD Part Number	Description	Size	Lead Length	Qty
SDRL-SE0045	Single, Pre-gelled, Green	45 x 35mm	1.0m	40 / Box
SDRL-SE0046	Single, Pre-gelled, Green	45 x 35mm	1.5m	40 / Box
SDRL-SE0047	Single, Pre-gelled, Green	45 x 35mm	2.5m	40 / Box

ARTG ID: 290956 | WAND: I70623-WAND-6OE2PQ



Surface Electrodes, Single Stimulating Sticky Pads, Disposable

SD Part Number	Description	Size	Lead Length	Qty
SDRL-SE0042	Small, Pre-gelled, Red, Blue, Black, Green	20 x 15mm	1.5m	80 / Box
SDRL-SE0043	Small, Pre-gelled, Red, Blue, Black, Green	20 x 15mm	2.5m	80 / Box
SDRL-SE0044	Small, Pre-gelled, Red, Blue, Black, Green	20 x 15mm	3.0m	80 / Box
SDRL-SE0040	Large, Pre-gelled, Red, Blue, Black, Green	20 x 25mm	1.5m	80 / Box
SDRL-SE0049	Large, Pre-gelled, Red, Blue, Black, Green	20 x 25mm	2.5m	80 / Box

ARTG ID: 290957 | WAND: I70623-WAND-6OE2PQ



Surface Electrodes, Parallel Pair Stimulating Sticky Pads, Disposable

SD Part Number	Description	Size	Lead Length	Qty
SDRL-SE0060	Small, Pre-gelled, Black / Red, Blue / Green	20 x 15mm	3.0m	20 Pair / Box
SDRL-SE0061	Small, Pre-gelled, Black / Red, Blue / Green	20 x 15mm	1.5m	20 Pair / Box
SDRL-SE0062	Small, Pre-gelled, Black / Red, Blue / Green	20 x 15mm	2.5m	40 Pair / Box
SDRL-SE0063	Large, Pre-gelled, Black / Red, Blue / Green	20 x 25mm	1.5m	40 Pair / Box
SDRL-SE1064	Large, Pre-gelled, Black / Red, Blue / Green	20 x 25mm	2.5m	40 Pair / Box

ARTG ID: 290957 | WAND: I70623-WAND-6OE2PQ



Surface Electrodes, Disposable



Surface Electrodes, Recording Sticky Pads, Single, Disposable

SD Part Number	Description	Size	Lead Length	Qty
SDRL-SE0050	Small, Pre-gelled, Red, Blue, Black, Green	20 x 15mm	1.0m	80 / Box
SDRL-SE0051	Small, Pre-gelled, Red, Blue, Black, Green	20 x 15mm	1.5m	80 / Box
SDRL-SE0052	Small, Pre-gelled, Red, Blue, Black, Green	20 x 15mm	2.5m	80 / Box
SDRL-SE0053	Large, Pre-gelled, Red, Blue, Black, Green	20 x 25mm	1.5m	80 / Box
SDRL-SE0054	Large, Pre-gelled, Red, Blue, Black, Green	20 x 25mm	2.5m	80 / Box

ARTG ID: 290957 | WAND: 170623-WAND-6OE2PQ

Surface Electrodes, Recording Sticky Pads, Parallel Pair, Disposable

SD Part Number	Description	Size	Lead Length	Qty
SDRL-SE0056	Small, Pre-gelled, Black / Red, Blue / Green	20 x 15mm	1.5m	20 Pair / Box
SDRL-SE0057	Small, Pre-gelled, Black / Red, Blue / Green	20 x 15mm	2.5m	20 Pair / Box
SDRL-SE0058	Small, Pre-gelled, Black / Red, Blue / Green	20 x 25mm	1.5m	40 Pair / Box
SDRL-SE0059	Large, Pre-gelled, Black / Red, Blue / Green	20 x 25mm	2.5m	40 Pair / Box

ARTG ID: 290957 | WAND: 170623-WAND-6OE2PQ

Extension Cables

SD Part Number	Description	Length	Qty
SDC-CC0023	Cadwell Dual Extension Lead Wire, 60" *	152cm	Each
SDC-CC0024	Cadwell Dual Extension Lead Wire, 180" *	457cm	Each
SDC-CC0025	Cadwell Quad Extension Lead Wire, 72" *	183cm	Each
SDC-CC0030	Reference Jumper. Two stackable 1.5mm, ultra-flexible *	15cm	Each
SDC-CC0029	Reference Jumper. Two stackable 1.5mm, ultra-flexible *	249cm	Each
SDFM-CC0019	Extension Lead Wire 6 Colours **	100cm	6 / Pk
SDFM-CC0020	Extension Lead Wire 6 Colours **	100cm	6 / Pk
SDFM-CC0021	Extension Lead Wire 6 Colours **	150cm	6 / Pk
SDFM-CC0022	Extension Lead Wire 6 Colours **	150cm	6 / Pk
SDFM-CC0023	Extension Lead Wire 6 Colours **	200cm	6 / Pk
SDFM-CC0024	Extension Lead Wire 6 Colours **	200cm	6 / Pk
SDFM-CC0025	Extension Lead Wire 6 Colours **	250cm	6 / Pk
SDFM-CC0026	Extension Lead Wire 6 Colours **	250cm	6 / Pk



Dual Extension Wire



SDC-CC0030




Generic Extension Lead Wire

* ARTG ID: 356742

** ARTG ID: 310991

CORTAC® Subdural Electrodes

PMT CORPORATION



Features

- Two contact material choices
 - Platinum (2110 series)
 - Stainless steel (2111 series)
- Multiple number of contacts (others available by non-standard request)
 - 2, 3, 4, 5, 6, 8, 10, 12, 14, 16, 18, 20, 24, 32, 36, 40, 48, 64
- Three standard contact sizes
 - 2.0 mm, 3.0 mm, 4.5 mm (1.0 mm available by non-standard request)
- 30+ grid configurations (columns x rows)
 - 2x2, 2x3, 2x4, 2x5, 2x6, 2x8, 3x2, 3x3, 3x4, 3x5, 3x6, 3x8, 4x2, 4x3, 4x4, 4x5, 4x6, 4x8, 5x2, 5x3, 5x4, 5x5, 5x6, 5x8, 6x2, 6x3, 6x4, 6x5, 6x6, 6x8, 8x2, 8x3, 8x4, 8x5, 8x6, 8x8
- Radiopaque option · Standard location is on contact #1
- Numbering format options · No numbers, top & bottom rows only, all contacts
- Contact spacing options (center to center) · 5 mm, 10 mm or 15 mm
- Three electrode thickness options on most electrodes · .5 mm, .76 mm, 1.01 mm
- Two surface options · smooth or enhanced (tx)
- Connector type options · in-line, mini (in-line is standard)
- Optional number of exits/tails · 1, 2, 3, 4, 6, 8
- Other options
 - Leuthardt grid perforated holes
 - Contour
 - L & t shaped strips and grids
 - Other options available upon request
 - Dual-sided/inter-hemispherical
 - Perpendicular tail design
 - Temporal lobe design

CORTAC® Cortical Electrodes, Strip (Sterile)

2mm, 3mm, or 4.5mm contact cortical electrodes based on total number of macro contacts.



SD Part Number	Description	Material	Qty	WAND
SDPM-SE0001	2 Contact Strip Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0018	2 Contact Strip Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ



SD Part Number	Description	Material	Qty	WAND
SDPM-SE0002	4 Contact Strip Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0019	4 Contact Strip Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ



SD Part Number	Description	Material	Qty	WAND
SDPM-SE0003	6 Contact Strip Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0020	6 Contact Strip Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ



SD Part Number	Description	Material	Qty	WAND
SDPM-SE0004	8 Contact Strip Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0021	8 Contact Strip Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ



SD Part Number	Description	Material	Qty	WAND
SDPM-SE0005	10 Contact Strip Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0022	10 Contact Strip Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ

CORTAC® Cortical Electrodes, Grid (Sterile)



SD Part Number	Description	Material	Qty	WAND
SDPM-SE0006	12 Contact Grid Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0023	12 Contact Grid Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ
SDPM-SE0007	16 Contact Grid Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0024	16 Contact Grid Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ
SDPM-SE0008	20 Contact Grid Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0025	20 Contact Grid Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ
SDPM-SE0009	22 Contact Grid (circle) Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0026	22 Contact Grid (circle) Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ
SDPM-SE0010	24 Contact Grid Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0027	24 Contact Grid Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ
SDPM-SE0011	32 Contact Grid Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0028	32 Contact Grid Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ
SDPM-SE0012	36 Contact Grid Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0029	36 Contact Grid Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ
SDPM-SE0013	40 Contact Grid Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0030	40 Contact Grid Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ
SDPM-SE0014	48 Contact Grid Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0031	48 Contact Grid Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ
SDPM-SE0015	64 Contact Grid Electrode	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0032	64 Contact Grid Electrode	Stainless Steel	Each	210119-WAND-6W6W60EJ
SDPM-SE0016	128 Contact Grid Electrode (assumed 1mm Contacts)	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0033	128 Contact Grid Electrode (assumed 1mm Contacts)	Stainless Steel	Each	210119-WAND-6W6W60EJ
SDPM-SE0017	256 Contact Grid Electrode (assumed 2mm Contacts)	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0034	256 Contact Grid Electrode (assumed 1mm Contacts)	Stainless Steel	Each	210119-WAND-6W6W60EJ



CORTAC® Subdural Electrodes, Strip (Sterile)



IOM One-Piece Design Platinum & Stainless Steel Electrodes

Supplied standard with 4.5mm contacts, but additional contact sizes and configurations available. Intra-Operative Monitoring Electrodes best for cases where the patient does not leave the OR.

SD Part Number	Description	Material	Qty	WAND
SDPM-SE0035	2 Contact Strip Electrode 1x2, 4.5mm Contacts, 6.5 foot tail with plugs	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0043	2 Contact Strip Electrode 1x2, 4.5mm Contacts, 6.5 feet Cable with 1.5mm plugs	Stainless Steel	Each	210119-WAND-6W6W60EJ




SD Part Number	Description	Material	Qty	WAND
SDPM-SE0036	4 Contact Strip Electrode 1x4, 4.5mm Contacts, 6.5 foot tail with plugs	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0044	4 Contact Strip Electrode 1x4, 4.5mm Contacts, 6.5 feet Cable with 1.5mm plugs	Stainless Steel	Each	210119-WAND-6W6W60EJ
SDPM-SE0037	4 Contact Strip Electrode 4x1, 4.5mm Contacts, 6.5 foot tail with plugs	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0045	4 Contact Strip Electrode 4x1, 4.5mm Contacts, 6.5 foot cable with 1.5mm plugs	Stainless Steel	Each	210119-WAND-6W6W60EJ



SD Part Number	Description	Material	Qty	WAND
SDPM-SE0038	6 Contact Strip Electrode 1x6, 4.5mm Contacts, 6.5 foot tail with plugs	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0046	6 Contact Strip Electrode 1x6, 4.5mm Contacts, 6.5 feet Cable with 1.5mm plugs	Stainless Steel	Each	210119-WAND-6W6W60EJ



SD Part Number	Description	Material	Qty	WAND
SDPM-SE0039	8 Contact Strip Electrode 1x8, 4.5mm Contact, 6.5 foot tail with plugs	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0047	8 Contact Strip Electrode 1x8, 4.5mm Contact, 6.5 feet Cable with 1.5mm plugs	Stainless Steel	Each	210119-WAND-6W6W60EJ



CORTAC® Subdural Electrodes, Grid (Sterile)

PMT® CORPORATION

IOM One-Piece Design Platinum & Stainless Steel Electrodes



SD Part Number	Description	Material	Qty	WAND
SDPM-SE0040	8 Contact Grid Electrode 2x4, 4.5mm Contacts, 6.5 foot tail with plugs	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0048	8 Contact Grid Electrode 2x4, 4.5mm Contacts, 6.5 feet Cable with 1.5mm plugs	Stainless Steel	Each	210119-WAND-6W6W60EJ




SD Part Number	Description	Material	Qty	WAND
SDPM-SE0041	12 Contact Grid Electrode 2x6, 4.5mm Contact, 6.5 foot tail with plugs	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0049	12 Contact Grid Electrode 2x6, 4.5mm Contact, 6.5 feet Cable with 1.5mm plugs	Stainless Steel	Each	210119-WAND-6W6W60EJ



SD Part Number	Description	Material	Qty	WAND
SDPM-SE0042	16 Contact Grid Electrode 4x4, 4.5mm Contacts, 6.5 foot tail with plugs	Platinum	Each	210119-WAND-6W6W60E4
SDPM-SE0050	16 Contact Grid Electrode 4x4, 4.5mm Contacts, 6.5 feet Cable with 1.5mm plugs	Stainless Steel	Each	210119-WAND-6W6W60EJ

Microwire Cortical Electrodes


Leuthardt electrodes have .76 mm grid thickness, 3.0 mm electrode disc diameter and perforated knock out holes for insertion of the microwire grommet. Other configurations available. Microwire depth electrodes - several configurations available as these electrodes are designed per customer spec.



SD Part Number	Description	Qty	WAND
SDPM-SC300	Macro-Micro Depth Electrode with 8 macros and 16 micros, several configurations available	Each	210119-WAND-6W60IP
SDPM-SC301	Macro-Micro Depth Electrode with 10 macros and 14 micros, several configurations available	Each	210119-WAND-6W60IP
SDPM-SC302	MicroBurst Microwire bundle with 10 microwires	Each	210119-WAND-6W60IP
SDPM-SC303	MicroBurst Macro Electrode with 4 macros	Each	210119-WAND-6W60IP
SDPM-SC304	MicroBurst Macro Electrode with 6 macros	Each	210119-WAND-6W60IP
SDPM-SC305	MicroBurst Macro Electrode with 8 macros	Each	210119-WAND-6W60IP
SDPM-SC306	MicroBurst Macro Electrode with 10 macros	Each	210119-WAND-6W60IP

Microwire Grommets





SD Part Number	Description	Material	Qty	WAND
SDPM-SE0150	16 contact (12 Recording, 4 Reference/ Ground) Microwire Array Electrode	Platinum	Each	210119-WAND-6W60IP
SDPM-SE0151	16 contact .75mm diameter Microwire Array Electrode	Platinum	Each	210119-WAND-6W60IP
SDPM-SE0152	New Non-standard Microwire Electrode	Platinum	Each	210119-WAND-6W60IP



sEEG Depth Electrodes

for 2.4mm Users (Sterile)

Mini Style Cables (2120) will be required in order to use sEEG Depth Electrodes. Stainless Steel stylet is provided. Contact "1" is located directly on the distal end of the electrode. Electrode contacts are available in 5mm sizing with 10mm spacing or 2mm contacts with 5mm Spacing.

SD Part Number	Description	Stylet	Material	Qty
SDPM-SE0053	4 Contact 12.5mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0054	6 Contact 19.5mm sEEG depth Electrode	Removable	Platinum	Each
				
SDPM-SE0055	8 Contact 26.5mm sEEG depth Electrode	Removable	Platinum	Each
				
SDPM-SE0056	10 Contact 33.5mm sEEG depth Electrode	Removable	Platinum	Each
				
SDPM-SE0057	12 Contact 40.5mm sEEG depth Electrode	Removable	Platinum	Each
				
SDPM-SE0058	14 Contact 47.5mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0059	16 Contact 54.5mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0060	16 Contact 61.5mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0061	16 Contact 68.5mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0062	16 Contact 80 mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0063	16 Contact 90 mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0064	4 Contact 17mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0065	6 Contact 27mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0066	8 Contact 37mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0067	10 Contact 47mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0068	12 Contact 57mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0069	14 contacts 67mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0070	16 Contact 77mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0071	16 Contact 34.3mm sEEG depth Electrode	Removable	Platinum	Each

WAND: 210119-WAND-6W60FS

sEEG Depth Electrodes

with Platinum Insulated Wiring (Sterile)

Available in 3.5mm CTC spacing. Mini Style Cables (2120) will be required in order to use sEEG Depth Electrodes.

SD Part Number	Description	Stylet	Material	Qty
SDPM-SE0072	4 Contact, 12.5mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0073	6 Contact, 19.5mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0074	8 Contact, 26.5mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0075	10 Contact, 33.5mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0076	12 Contact, 40.5mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0077	14 Contact, 47.5mm sEEG depth Electrode	Removable	Platinum	Each
SDPM-SE0078	16 Contact, 54.5mm sEEG depth Electrode	Removable	Platinum	Each

WAND: 210119-WAND-6W60FS

DEPTHALON Depth Electrodes

1.2mm diameter and more rigid, can be used with or without stereotactic placement.

Platinum (Sterile)

SD Part Number	Description	Diameter	Material	Qty
SDPM-SE0079	4 Contact Platinum Depth Electrode	1.2mm	Platinum	Each
SDPM-SE0080	6 Contact Platinum Depth Electrode	1.2mm	Platinum	Each
SDPM-SE0081	8 Contact Platinum Depth Electrode	1.2mm	Platinum	Each
SDPM-SE0082	10 Contact Platinum Depth Electrode	1.2mm	Platinum	Each

WAND: 210119-WAND-6W60FS

Stainless Steel (Sterile)

SD Part Number	Description	Diameter	Material	Qty
SDPM-SE0083	4 Contact Stainless Steel Depth Electrode	1.2mm	Stainless Steel	Each
SDPM-SE0084	6 Contact Stainless Steel Depth Electrode	1.2mm	Stainless Steel	Each
SDPM-SE0085	8 Contact Stainless Steel Depth Electrode	1.2mm	Stainless Steel	Each
SDPM-SE0086	10 Contact Stainless Steel Depth Electrode	1.2mm	Stainless Steel	Each

WAND: 210119-WAND-6W60G5

ACCUNAV[®] Neuronavigation Depth Electrode, (Sterile)

PMT[®] CORPORATION

PMT ACCUNAV[®] depth electrodes have been developed to allow surgeons the ultimate advancement in epilepsy electrode intracranial placement. The ACCUNAV[®] is an all new depth electrode that was specifically designed and produced for use with electromagnetic and optical neuronavigation systems.

Features

- Internal lumen provides ease of neuronavigation stylet placement within the electrode.
 - Smooth outer tubing design with integrated electrodes for reduced brain trauma.
 - Extra malleable substrate tubing material reduces the potential for hemotoma during placement and duration of monitoring.
 - ACCUNAV[®] is 10% smaller in diameter than all other neuronavigation depth electrodes on the market.
 - Configurations are available with 2, 4, 6 or 8 contacts, with two spacing options providing 50% more recording depth than other neuronavigation depths offered on the market.
 - Distal ends records to the very tip
 - Electrode tails connect to existing PMT inline connectors.
- PMT ACCUNAV[®] depth electrodes are available with 2, 4, 6 or 8 contacts.
 - All models have an internal lumen designed to work with the two major neuronavigation systems currently on the market in the us.
 - All models are 1.8mm in diameter and work with any stereotactic system's guiding stylet, accommodates a Medtronic Axiom Stylet, or Brainlab optical stylet.
 - All contacts are 2mm in length with 5mm or 10mm center to center spacing.



SD Part Number	Description	Qty
SDPM-SE0087	2 Contacts, 7mm Recording Depth, 5mm Spacing	Each
SDPM-SE0088	2 Contacts, 12mm Recording Depth, 10mm Spacing	Each
SDPM-SE0089	4 Contacts, 17mm Recording Depth, 5mm Spacing	Each
SDPM-SE0090	4 Contacts, 32mm Recording Depth, 10mm Spacing	Each
SDPM-SE0091	6 Contacts, 27mm Recording Depth, 5mm Spacing	Each
SDPM-SE0092	6 Contacts, 52mm Recording Depth, 10mm Spacing	Each
SDPM-SE0093	8 Contacts, 37mm Recording Depth, 5mm Spacing	Each
SDPM-SE0094	8 Contacts, 72mm Recording Depth, 10mm Spacing	Each

WAND: 210119-WAND-6W601P

Cables, Kits & Spares

Interconnection Cables

SD Part Number	Description	Qty
SDPM-CC0001	Internconnection Cable, Mini Connector, 4 Contact	Each
SDPM-CC0002	Internconnection Cable, Mini Connector, 6 Contact	Each
SDPM-CC0003	Internconnection Cable, Mini Connector, 8 Contact	Each
SDPM-CC0004	Internconnection Cable, Mini Connector, 10 Contact	Each
SDPM-CC0005	Internconnection Cable, Mini Connector, 12 Contact	Each
SDPM-CC0006	Internconnection Cable, Mini Connector, 14 Contact	Each
SDPM-CC0007	Internconnection Cable, Mini Connector, 16 Contact	Each
SDPM-CC0008	Modular Testing Cable Box for Mini style tails, 4-16 contacts	Each
SDPM-CC0009	Internconnection Cable, Inline Connector, 4 contact	Each
SDPM-CC0010	Internconnection Cable, Inline Connector, 8 contact	Each
SDPM-CC0011	Internconnection Cable, Inline Connector, 12 contact	Each
SDPM-CC0012	Internconnection Cable, Inline Connector, 16 contact	Each
SDPM-CC0013	Internconnection Cable, Inline Connector, 6 contact	Each
SDPM-CC0014	Internconnection Cable, Inline Connector, 12 contact	Each
SDPM-CC0015	Internconnection Cable, Inline Connector, 24 contact	Each
SDPM-CC0016	Internconnection Cable, Inline Connector, 8 contact	Each
SDPM-CC0017	Internconnection Cable, Inline Connector, 16 contact	Each
SDPM-CC0018	Internconnection Cable, Inline Connector, 24 contact	Each
SDPM-CC0019	Internconnection Cable, Inline Connector, 32 contact	Each
SDPM-CC0020	Internconnection Cable, Inline Connector, 48 contact	Each
SDPM-CC0021	Internconnection Cable, Inline Connector, 64 contact	Each
SDPM-CC0022	Internconnection Cable, Inline Connector, 10 contact	Each
SDPM-CC0023	Internconnection Cable, Inline Connector, 20 contact	Each
SDPM-CC0024	Internconnection Cable, Inline Connector, 40 contact	Each
SDPM-CC0025	Internconnection Cable, Inline Connector, 12 contact	Each
SDPM-CC0026	Internconnection Cable, Inline Connector, 24 contact	Each
SDPM-CC0027	Internconnection Cable, Inline Connector, 36 contact	Each
SDPM-CC0028	Internconnection Cable, Inline Connector, 48 contact	Each
SDPM-CC0029	Internconnection Cable, Inline Connector, 14 contact	Each
SDPM-CC0030	Internconnection Cable, Inline Connector, 16 contact	Each
SDPM-CC0031	Internconnection Cable, Inline Connector, 32 contact	Each
SDPM-CC0032	Internconnection Cable, Inline Connector, 48 contact	Each
SDPM-CC0033	Internconnection Cable, Inline Connector, 64 contact	Each



Additional configurations, like Quick Disconnect are available

WAND: 210119-WAND-6W60HI

Cables, Kits & Spares

Frame Specific Kits

Additional frame specific kits available (FHC, Neuromate, etc.).

SD Part Number	Description	Qty
SDPM-K0001	2.7mm Leksell Kit: Includes 2103-26, 39, 58, 59, 62, 61, 73, 74	Kit
SDPM-K0002	3.2mm Leksell Kit: Includes 2103-26, 39, 58, 59, 62, 61, 76, 77	Kit
SDPM-K0003	2.1mm Leksell Kit: Includes 2103-142, 52, 92, 60, 68, 69, 87, 88, 93, 97	Kit
SDPM-K0004	2.4mm Leksell Kit: Includes 2103-143, 52, 92, 63, 65, 66, 69, 107, 108	Kit
SDPM-K0005	2.4mm ROSA Kit: Includes 2103-52, 63, 65, 66, 69, 92, 143 (2)	Kit
SDPM-K0006	2.1mm ROSA Kit: Includes 2103-142 (2), 52, 60, 68, 69, 85, 92, 93, 94	Kit
SDPM-K0007	2.4mm Brainlab Varioguide Kit: Includes 2103-42, 143 (2), 52, 92, 63, 66, 69, 65	Kit
SDPM-K0008	2.1mm Brainlab Varioguide Kit: Includes 2103-71, 142 (2), 52, 92, 60, 68, 69, 93	Kit
SDPM-K0009	2.7mm Stealth Vertek kit: Includes 2103-56, 102, 103	Kit
SDPM-K0010	3.2mm Stealth Vertek kit: Includes 2103-80, 102, 103	Kit
SDPM-K0011	2.1 Stealth Vertek kit: Includes 2103-142, 52, 60, 03, 68, 69, 85, 92, 93	Kit
SDPM-K0012	2.4 Stealth Vertek kit: Includes 2103-143, 52, 63, 66, 69, 92, 65, 114	Kit
SDPM-K0013	2.7mm BRW/CRW Kit: Includes 2103-26, 40, 45, 47, 48, 51	Kit
SDPM-K0014	3.2mm BRW/CRW Kit: Includes 2103-26, 33, 40, 45, 47, 51	Kit
SDPM-K0015	2.1mm BRW/CRW Kit: includes 2103-17, 34, 52, 92, 60, 68, 69, 85, 93	Kit
SDPM-K0016	2.4mm BRW/CRW Kit: Includes 2103-36, 95, 52, 92, 63, 66, 69, 65	Kit

WAND: 210119-WAND-6W60IP

Cables, Kits & Spares

Frame Specific Guide Bushings - sEEG

Additional bushings available (FHC, Neuromate, other non-standard designs).

SD Part Number	Description	Qty
SDPM-MI0001	Leksell upper guide bushing for 2.1mm drill	Each
SDPM-MI0002	Leksell lower guide bushing for 2.1mm drill	Each
SDPM-MI0003	Leksell-sEEG drill reducing tube for 2.1mm	Each
SDPM-MI0004	Leksell upper guide bushing for 2.4mm drill	Each
SDPM-MI0005	Leksell lower guide bushing for 2.4mm drill	Each
SDPM-MI0006	ROSA guide bushing-long for 2.1mm drill (standard)	Each
SDPM-MI0007	ROSA guide reducing tube-long for 2.1mm drill	Each
SDPM-MI0008	ROSA guide bushing-short for 2.1mm drill (non-standard)	Each
SDPM-MI0009	ROSA guide reducing tube- short for 2.1mm drill	Each
SDPM-MI0010	ROSA guide bushing- long for 2.4mm drill (standard)	Each
SDPM-MI0011	ROSA guide bushing-short for 2.4mm drill (non-standard)	Each
SDPM-MI0012	Brainlab Varioguide guide bushing for 2.1mm drill	Each
SDPM-MI0013	Brainlab Varioguide reducing tube for 2.1mm drill	Each
SDPM-MI0014	Brianlab Varioguide guide bushing for 2.4mm drill	Each
SDPM-MI0015	Stealth Vertek guide bushing for 2.1mm drill	Each
SDPM-MI0016	Stealth Vertek guide reducing tube for 2.1mm drill	Each
SDPM-MI0017	Stealth Vertek guide bushing for 2.4mm drill	Each
SDPM-MI0018	BRW/ CRW guide bushing for 2.1mm drill	Each
SDPM-MI0019	BRW/CRW reducing tube for 2.1mm drill	Each
SDPM-MI0020	BRW/CRW guide bushing for 2.4mm drill	Each

WAND: 210119-WAND-6W60IP

Cables, Kits & Spares

Universal sEEG Kit Parts

SD Part Number	Description	Qty
SDPM-MI0021	Depth Ruler for sEEG depth	Each
SDPM-MI0022	Guiding Stylet for sEEG depth placement	Each
SDPM-MI0023	Depth stop for sEEG guiding stylet	Each
SDPM-MI0024	sEEG depth anchor bolt driver for 2.1mm	Each
SDPM-MI0025	sEEG anchor bolt removal driver 1" length for 2.1mm	Each
SDPM-MI0026	Depth Stop for 2.1mm drill bit for sEEG anchor bolt	Each
SDPM-MI0027	Depth stop for sEEG anchor bolt driver for 2.1mm	Each
SDPM-MI0028	Depth stop for 2.4mm drill bit for sEEG anchor bolt for 2.4mm	Each
SDPM-MI0029	Depth stop for sEEG anchor bolt driver for 2.4mm	Each
SDPM-MI0030	sEEG depth anchor bolt driver for 2.4mm	Each
SDPM-MI0031	sEEG anchor bolt removal driver 1" length for 2.4mm	Each
SDPM-MI0032	Circular Depth Stop for 2.4mm Drill Bit	Each
SDPM-MI0033	Allen Wrench for Circular Depth Stop for 2.4mm Drill Bit	Each

WAND: 210119-WAND-6W60IP

Universal Standard Depth and Visualase Kit Parts

SD Part Number	Description	Qty
SDPM-MI0060	T-Driver	Each
SDPM-MI0061	Depth Stop for 3.2mm Drill Bit	Each
SDPM-MI0062	Depth Stop for 2103-26 Anchor Bolt Driver	Each
SDPM-MI0063	Guiding stylet, 1.2mm diameter, for depth and microwire depth electrodes	Each
SDPM-MI0064	Depth Stop for 1.2mm stylet	Each

WAND: 210119-WAND-6W60IP

Cables, Kits & Spares

Frame Specific Guide Bushing - Standard Depth and Visualase

SD Part Number	Description	Qty
SDPM-MI0034	Leksell lower guide bushing for 2.7mm drill	Each
SDPM-MI0035	Leksell upper guide bushing for 2.7mm drill	Each
SDPM-MI0036	Leksell lower guide bushing for 3.2mm drill	Each
SDPM-MI0037	Leksell upper guide bushing for 3.2mm drill	Each
SDPM-MI0038	Leksell Upper Split Cannula Bushing for T-Driver	Each
SDPM-MI0039	Leksell Lower Split Cannula Bushing for T-Driver	Each
SDPM-MI0040	Leksell Upper Split Cannula Bushing	Each
SDPM-MI0041	Leksell Lower Split Cannula Bushing	Each
SDPM-MI0042	Leksell Split Cannula	Each
SDPM-MI0043	Stealth/Vertek Drill bushing for 2.7mm Drill	Each
SDPM-MI0044	Stealth/Vertek Drill bushing for 3.2mm Drill	Each
SDPM-MI0045	Stealth/Vertek Split Bushing for Standard Depth Anchor Bolt Placement	Each
SDPM-MI0046	ROSA Guide Bushing for depthalon T-Driver for either 2.7 or 3.2mm drills	Each
SDPM-MI0047	ROSA Reducing bushing w/ thumbscrew for 3.2mm drill	Each
SDPM-MI0048	ROSA Reducing bushing w/ thumbscrew for 2.7mm drill	Each
SDPM-MI0049	Brainlab VARIOguide 2 piece SPLIT bushing	Each
SDPM-MI0050	Brainlab VARIOguide ID 4.0mm OD drill guide for 2.7mm drill	Each
SDPM-MI0051	Brainlab VARIOguide ID 4.0mm OD drill guide for 3.2mm drill	Each
SDPM-MI0052	Brainlab VARIOguide 2 piece bushing for depth	Each
SDPM-MI0053	BRW/CRW Split Cannula length 240mm	Each
SDPM-MI0054	BRW/CRW Guide Bushing for T-Driver	Each
SDPM-MI0055	Drill guide bushing for 3.2mm drill	Each
SDPM-MI0056	BRW/CRW Split cannula bushing	Each
SDPM-MI0057	BRW/CRW guide bushing for 2.7mm drill	Each
SDPM-MI0058	Thumb screw (For 2103-2)	Each
SDPM-MI0059	BRW/CRW guide bushing for 3.2mm drill	Each

WAND: 210119-WAND-6W60IP

Cables, Kits & Spares

Misc. Components for Depth and sEEG


SD Part Number	Description	Qty
SDPM-MI0065	PMT Anchor bolt for 3.2mm Standard depth (available in 25mm, 30mm, 35mm, and 40mm)	Each
SDPM-MI0066	PMT Anchor bolt for 2.7mm Standard depth (available in 25mm, 30mm, 35mm, and 40mm)	Each
SDPM-MI0067	PMT Anchor Bolt for 2.7mm drill and laser probe compatible (available in 25mm, 30mm, 35mm, and 40mm)	Each
SDPM-MI0068	PMT Anchor Bolt for 3.2mm drill and laser probe compatible (available in 25mm, 30mm, 35mm, and 40mm)	Each
SDPM-MI0069	PMT Anchor Bolt Cap for Standard Depth Electrode (sold as extra, comes with anchor bolt)	Each
SDPM-MI0070	PMT Anchor bolt reducing system for laser probe compatible bolt	Each
SDPM-MI0071	PMT sEEG Anchor Bolt for 2.1mm Drill (available in 20mm, 25mm, 30mm, 35mm, and 40mm)	Each
SDPM-MI0072	PMT sEEG Anchor Bolt Cap for 2.1mm Drill (2103-15-XX Anchor Bolt) sold as extra, comes with sEEG	Each
SDPM-MI0073	PMT sEEG Anchor Bolt for 2.4mm Drill (available in 20mm, 25mm, 30mm, 35mm, and 40mm)	Each
SDPM-MI0074	PMT sEEG Anchor Bolt Cap for 2.4mm Drill (2103-29-XX Anchor Bolt) sold as extra, comes with sEEG	Each
SDPM-MI0075	PMT sEEG Anchor Bolt Cap for 2.4mm Drill (2103-29-XX Anchor Bolt)	Each
SDPM-MI0076	PMT sEEG Anchor Bolt Cap Cover- Yellow- for folding over tails	Each
SDPM-MI0077	PMT sEEG Rescue Anchor Bolt for 2.4mm drill (available in 25mm, 30mm, and 35mm)	Each
SDPM-MI0078	PMT sEEG Rescue Anchor Bolt for 2.1mm drill (available in 25mm, 30mm, and 35mm)	Each
SDPM-MI0079	Twist Drill Bit 2.4mm x 254mm (10in.) Sterile, Single use	Each
SDPM-MI0080	Twist Drill Bit 2.1mm x 254mm (10in.) Sterile, Single use	Each
SDPM-MI0081	In-line needles, 13GA, (no charge with in-line tail electrode order)	Each
SDPM-MI0082	CRW Sterilization Tray for sEEG 2.1mm bit	Each
SDPM-MI0083	CRW Sterilization Tray for sEEG 2.4mm bit	Each
SDPM-MI0084	Brainlab Sterilization Tray for sEEG 2.1mm bit	Each
SDPM-MI0085	Brainlab Sterilization Tray for sEEG 2.4mm bit	Each
SDPM-MI0086	Leksell Sterilization Tray for sEEG 2.1mm bit	Each
SDPM-MI0087	Leksell Sterilization Tray for sEEG 2.4mm bit	Each
SDPM-MI0088	Rosa Sterilization Tray for sEEG 2.1mm bit	Each
SDPM-MI0089	Rosa Sterilization Tray for sEEG 2.4mm bit	Each
SDPM-MI0090	Stealth Sterilization Tray for sEEG 2.1mm bit	Each
SDPM-MI0091	Stealth Sterilization Tray for sEEG 2.4mm bit	Each

WAND: 210119-WAND-6W60IP

Sphenoidal Electrodes

Sphenoidal Electrodes Kits

Includes two (2) electrodes, two (2) needles and two (2) cables
 Can be supplied individually or in sets, using 36 Gauge Platinum, Stainless Steel or Silver wire unless specified otherwise sphenoidal electrodes arrive in 30 cm length, with 3-5mm tip.



SD Part Number	Description	Qty
SDPM-K0017	2 Each Stainless Steel Electrodes, 2 Needles and 2 Cables - Crimped	Kit
SDPM-K0018	2 Each Stainless Steel Electrodes, 2 Needles and 2 Cables - Non crimped	Kit
WAND: 210119-WAND-6W6W60EJ		
SDPM-K0019	2 Each Platinum Electrodes, 2 Needles and 2 Cables - Crimped	Kit
SDPM-K0020	2 Each Platinum Electrodes, 2 Needles and 2 Cables - Non crimped	Kit
WAND: 210119-WAND-6W6W60E4		

Sphenoidal Electrodes Sets

Include two (2) electrodes and two (2) needles

SD Part Number	Description	Qty
SDPM-K0021	2 Each Stainless Steel Electrodes and 2 Needles - Crimped	Kit
SDPM-K0022	2 Each Stainless Steel Electrodes and 2 Needles - Non crimped	Kit
WAND: 210119-WAND-6W6W60EJ		
SDPM-K0023	2 Each Platinum Electrodes and 2 Needles - Crimped	Kit
SDPM-K0024	2 Each Platinum Electrodes and 2 Needles - Non crimped	Kit
WAND: 210119-WAND-6W6W60E4		

Sphenoidal Electrode additional components

SD Part Number	Description	Qty
SDPM-CC0035	Sphenoidal Electrode EEG Cable 1.5mm socket Set (2 cables)	Each
SDPM-CC0036	Crimping Tool for Non-Crimped Sphenoidal Electrode	Each
WAND: 210119-WAND-6W6W60EJ		



Symbiotic Devices was founded in 2012 on a passion not only for neuroscience and neurotechnology, but also to provide clients with uncompromised quality, reliability and assistance. Striving to contribute to a sustainable successful research and clinical sectors across Australia and New Zealand, Symbiotic Devices aims to balance quality, accessibility and customisation.

We know that each neurophysiology project will present unique challenges, and our wide range of partners means we are able to provide customised and dynamic solutions to fit our client's needs.

First partnering with Brain Products, EasyCap and Rogue Research, Symbiotic Devices now proudly represents CED, Deymed, Cadwell, NeuroConn/NeuroCare, Rhythmlink, Pearl Technology and VPixx.

If you want the best solutions for your neuroscience and neurophysiology needs, you can trust Symbiotic Devices to provide the highest quality equipment, applications and support while contributing to a sustainable and ethical future for neuroscience practices.



SYMBIOTIC
DEVICES

au 1300 934 947 nz 09 886 4877
team@symbioticdevices.com.au

Unit 6, 105-111 Ricketts Rd,
Mount Waverley, VIC 3149

www.symbioticdevices.com.au



Contact us

team@symbioticdevices.com.au | Ph: 1300 934 947

SD-B0007-NZ